

Selected FastTag[®] references:

Graber TE, Baird SD, Kao PN, Mathews MB, and Holcik M. 2010. NF45 functions as an IRES *trans*-acting factor that is required for translation of cIAP1 during the unfolded protein response. *Cell Death and Differentiation*. 17:719-29.

Harris TJ, Green JJ, Fung PW, Langer R, Anderson DG, and Bhatia SN. 2010. Tissue-specific gene delivery via nanoparticle coating. *Biomaterials*. 31(5):998-1006.

Valentin D, Marin A, Hedge A-M, and Rowe PSN. 2009. Matrix extracellular phosphoglycoprotein (MEPE) is a new bone renal hormone and vascularization modulator. *Endocrinol*. 150:4012-4023.

Valentin D, Marin A, Hedge A-M, and Rowe PSN. 2009. MEPE is a new bone renal hormone and vascularization modulator. *Endocrinol*. Doi:10.1210/en.2009-0216.

Grosse S, Thevenot G, Aron Y, Duverger E, Abdelkarim M, Roche A-C, Monsigny M, Fajac I. 2008. *In vivo* gene delivery in the mouse lung with lactosylated polyethylenimine, questioning the relevance of *in vitro* experiments. *J Control. Release*. 132:105-12

Shiraishi Y, Asano K, Niimi K, Fukunaga K, Wakaki M, Kagyo J, Takihara T, Ueda S, Nakajima T, Oguma T, Suzuki Y, Shiomi T, Sayama K, Kagawa S, Ikeda E, Hirai H, Nagata K, Nakamura M, Miyasho T, Ishizaka A. 2008. Cyclooxygenase-2/Prostaglandin D2/CRTH2 pathway mediates double-stranded RNA-induced enhancement of allergic airway inflammation. *J Immunol*. 180:541-49.

Grosse S, Aron Y, Thevenot G, Monsigny M, Fajac I. 2007. Cytoskeletal involvement in the cellular trafficking of plasmid/PEI derivative complexes. *J Control. Release*. 122:111-17.

Niimi K, Asano K, Shiraishi Y, Nakajima T, Wakaki M, Kagyo J, Takihara T, Suzuki Y, Fukunaga K, Shiomi T, Oguma T, Sayama K, Yamaguchi K, Natori Y, Matsumoto M, Seya T, Yamaya M, Ishizaka A. 2007. TLR3-mediated synthesis and release of eotaxin-1/CCL11 from human bronchial smooth muscle cells stimulated with double-stranded RNA. *J Immunol*. 178:489-95.

Shimada S, Sunagawa M, Hanashiro K, Nakamura M, Kosugi T. 2007. RNA interference targeting embryonic myosin heavy chain isoform inhibited mRNA expressions of phenotype markers in rabbit cultured vascular smooth muscle cells. *Heart Vessels*. 22:41-47.

Fule T, Mathe M, Suba Z, Csapo Z, Szarvas T, Tatrai P, Paku S, Kovalszky I. 2006. The presence of human papillomavirus 16 in neural structures and vascular endothelial cells. *Virology*. 348:289-96.

Grosse S, Thevenot G, Monsigny M, Fajac I. 2006. Which mechanism for nuclear import of plasmid DNA complexed with polyethylenimine derivatives? *J Gene Medicine* 8(7):845-51.

Funhoff AM, van Nostrum CF, Lok MC, Kruijtz JAW, Crommelin DFA, Hennink WE. 2005. Cationic polymethacrylates with covalently linked membrane destabilizing peptides as gene delivery vectors. *J Controlled Release*. 101:233-46.

Gasiorowski JZ and Dean DA. 2005. Postmitotic nuclear retention of episomal plasmids is altered by DNA labeling and detection methods. *Mol. Therapy*. 12(3):460-67.

Grosse S, Aron Y, Thevenot G, Francois D, Monsigny M, Fajac I. 2005. Potocytosis and cellular exit of complexes as cellular pathways for gene delivery by polycations. *J Gene Med*. 7:1275-86.

Kawase A, Kobayashi N, Isaji K, Nishikawa M, Takakura Y. 2005. Manipulation of local disposition and gene expression characteristics of plasmid DNA following intramuscular administration by complexation with cationic macromolecule. *Int J Pharmaceutics*. 293:291-301.

Kwak S-E, Kim J-E, Kim D-S, Won M-H, Choi H-C, Kim Y-I, Song H-K, Choi S-Y, Kang T-C. 2005. Differential effects of vigabatrin and zonisamide on the neuropeptide Y system in the hippocampus of seizure prone gerbil. *Neuropeptides* 39(5):507-13.

Yamane I, Nishikawa M, Takakura Y. 2005. Cellular uptake and activation characteristics of naked plasmid DNA and its cationic liposome complex in human macrophages. *Int J Pharmaceutics*. 305(1-2):145-53.

Yasuda K, Ogawa Y, Yamane I, Nishikawa M, Takakura Y. 2005. Macrophage activation by a DNA/cationic liposome complex requires endosomal acidification and TLR9-dependent and -independent pathways. *J Leukocyte Biol*. 77:71-9.

Zheng P-P. 2005. Low-molecular weight caldesmon in glioma neovascularization. Ph.D. Thesis. <http://www.ep.eur.nl/bitstream/1765/6736/1/thesis+zheng+drukker+definitieve+versie+3+jan+2005.pdf>

Funhoff AM, van Nostrum CF, Koning GA, Schuurmans-Nieuwenbroek NM, Crommelin DJ, Hennink WE. 2004. Endosomal escape of polymeric gene delivery complexes is not always enhanced by polymers buffering at low pH. *Biomacromolecules*. 5(1):32-9.

Funhoff AM, van Nostrum CF, Lok MC, Fretz MM, Crommelin DFA, Hennink WE. 2004. Poly(3-guanidinopropyl methacrylate): A novel cationic polymer for gene delivery. *Bioconjugate Chem*. 15:1212-20.

Grosse S, Aron Y, Honore I, Thevonot G, Danel C, Roch AC, Monsigny M, Fajac I.

2004. Lactosylated polyethylenimine for gene transfer into airway epithelial cells: Role of the sugar moiety in cell delivery and intracellular trafficking of the complexes. *J Gene Med.* 6(3):345-56.

Hisazumi J, Kobayashi N, Nishikawa M, Takakura Y. 2004. Significant role of liver sinusoidal endothelial cells in hepatic uptake and degradation of naked plasmid DNA after intravenous infection. *Pharm Res.* 21(7):1223-8.

Wiethoff CM, Koe JG, Koe GS, Middaugh CR. 2004. Compositional effects of cationic lipid/DNA delivery systems on transgene expression in cell culture. *J Pharm Sci.* 93(1):108-23.

Yasuda K, Kawano H, Yamane I, Ogawa Y, Yoshinaga T, Nishikawa M, Takakura Y. 2004. Restricted cytokine production from mouse peritoneal macrophages in culture in spite of extensive uptake of plasmid DNA. *Immunology.* 111:282-90.

Rejman J, Wagenaar A, Engberts J, Hoekstra D. 2004. Characterization and transfection properties of lipoplexes stabilized with novel exchangeable polyethylene glycol-lipid conjugates. *Biochem. Biophys. Acta.* 1660:41-52.

Kawase A, Nomura T, Yasuda K, Kobayashi N, Hashida M, Takakura Y. 2003. Disposition and gene expression characteristics in solid tumors and skeletal muscle after direct injection of naked plasmid DNA in mice. *J Pharm Sci.* 92(6):1295-304.

Luten J, van Steenis JH, van Someren R, Kemmink J, Schuurmans-Nieuwenbroeck NM, Koning GA, Crommelin DJ, van Nostrum CF, Hennink WE. 2003. Water-soluble biodegradable cationic polyphosphazenes for gene delivery. *J Controlled Release.* 89(3):483-97.

Silvestrini G, Mocetti P, Di Grezia R, Berni S, Bonucci E. 2003. Localization of the glucocorticoid receptor mRNA in cartilage and bone cells of the rat. An *in situ* hybridization study. *Eur J Histochem.* 47(3):245-52.

Grosse S, Tremeau-Bravard A, Aron Y, Briand P, Fajac I. 2002. Intracellular rate-limiting steps of gene transfer using glycosylated polylysines in cystic fibrosis airway epithelial cells. *Gene Ther.* 9(15):1000-07.

Yoshinaga T, Yasuda K, Ogawa Y, Takakura Y. 2002. Efficient uptake and rapid degradation of plasmid DNA by murine dendritic cells via a specific mechanism. *J. Biochem. Biophys. Res. Comm.* 299:389-94.

Kobayashi N, Kuramoto T, Yamaoka K, Hashida M, Takakura Y. 2001. Hepatic uptake and gene expression mechanisms following intravenous administration of plasmid DNA by conventional and hydrodynamics-based procedures. *J. Pharm. Exp Ther.* 297(3):853-60.

Nishikawa M, Yamauchi M, Morimoto K, Ishida E, Takakura Y, Hashida M. 2000. Hepatocyte-targeted *in vivo* gene expression by intravenous injection of plasmid DNA complexed with synthetic multi-functional gene delivery system. *Gene Therapy*. 7:548-55.

Sakurai F, Inoue R, Nishino Y, Okuda A, Matsumoto O, Taga T, Yanashita F, Takakura Y, Hashida M. 2000. Effect of DNA/liposome mixing ratio on the physicochemical characteristics, cellular uptake and intracellular trafficking of plasmid DNA/cationic liposome complexes and subsequent gene expression. *J. Control Release*. 66:255-69.

Ely SA, Powers J, Lewis D, Chang S, Rubio A, O'Leary J., Knowles DM. 1999. Kaposi's sarcoma-associated herpesvirus-positive primary effusion lymphoma arising in the subarachnoid space. *Human Pathology*. 30(8):981-84.

MacDonald RC, Rakhmanova VA, Choi KL, Rosenzweig HS, Lahiri MK. 1999. O-ethylphosphatidylcholine: A metabolizable cationic phospholipid which is a serum-compatible DNA transfection agent. *J Pharm Sci*. 88(9):896-904.

Pitard B, Oudrhiri N, Vigneron J-P, Hauchecorne M, Aguerre A, Toury R, Airiau M, Ramasawmy R, Scherman D, Crouzet J, Lehn J-M, Lehn P. 1999. Structural characteristics of supramolecular assemblies formed by guanidinium-cholesterol reagents for gene transfection. *Proc. Nat. Acad. Sci*. 96(6):2621-26.

Takakura Y, Takagi T, Hashiguchi M, Nishikawa M, Yamashita F, Doi T, Imanishi T, Suzuki H, Kodama T, Hashida M. 1999. Characterization of Plasmid DNA Binding and Uptake by Peritoneal Macrophages from Class A Scavenger Receptor Knockout Mice. *Pharmaceut. Res*. 16:503-08.

Ace CI and Okulicz WC. 1998. A progesterone-induced endometrial homolog of a new candidate tumor suppressor, DMBT1. *J. Clin. Endocrin. and Metab*. 83(10):3569-73.

Daniel SG, Westling ME, Moss MS, Kanagy BD. 1998. FastTagTM nucleic acid labeling system: A versatile method for incorporating haptens, fluorochromes, and affinity ligands into DNA, RNA, and oligonucleotides. *BioTechniques*. 24:484-89.

Jung J, Moronne MM, and Weier H-UG. 1998. Simple blocking system for use with maleimide-labeled nucleic acid probes. *BioTechniques*. 25:390-94.

Reed JA, et al. 1998. Demonstration of Kaposi's sarcoma-associated herpes virus cyclin D homolog in cutaneous Kaposi's sarcoma by colorimetric *in situ* hybridization using a catalyzed signal amplification system. *Blood*. 91:3825-32.

Takagi T, et al. 1998. Involvement of specific mechanism in plasmid DNA uptake by mouse peritoneal macrophages. *Biochem Biophys Res Commun*. 245:729-33.